

KAB Dental Inc.

Ph : 586-983-2502, 800-422-3520 • FAX : 586 983 2506

HANDPIECE WORK ORDER

OFFICE NAME		PHONE	FAX	
ADDRESS		CITY	STATE	ZIP CODE
PAYMENT INFORMATION	<input type="checkbox"/> INVOICE / PAY ON RECEIPT	<input type="checkbox"/> VISA	<input type="checkbox"/> MASTERCARD	CARD NUMBER
				EXP. DATE

1	MAKE/MODEL	SERIAL NUMBER	REPAIR OPTION (SEE BELOW)
	COMMENTS (NATURE OF PROBLEM, DETAILS, ETC.)		
2	MAKE/MODEL	SERIAL NUMBER	REPAIR OPTION
	COMMENTS		
3	MAKE/MODEL	SERIAL NUMBER	REPAIR OPTION
	COMMENTS		
4	MAKE/MODEL	SERIAL NUMBER	REPAIR OPTION
	COMMENTS		

PLEASE RECORD ONE HANDPIECE PER LINE ABOVE. INCLUDE ADDITIONAL HANDPIECE INFORMATION ON A SECOND FORM OR PIECE OF PAPER.

DATE	REQUESTED BY	SIGNATURE
------	--------------	-----------

HANDPIECE REPAIR OPTIONS


HIGH SPEED HANDPIECES

OPT.	DESCRIPTION	
A***	STANDARD TURBINE OVERHAUL : disassembly of turbine, cleaning, and replacement of all wear parts (bearings, o-rings, etc); 4 month warranty	\$88.95
	CANISTER-STYLE HANDPIECES : replacement turbine canister assembly; 90 day warranty	\$66.95
	PUSHBUTTON CANISTER-STYLE HANDPIECES : replacement pushbutton-type turbine canister assembly; 90 day warranty	\$85.00
	STAR LUBE-FREE HANDPIECES : disassembly of turbine, cleaning, and replacement of all wear parts (bearings, o-rings, etc) using "maintenance-free" ceramic bearings; 6 month warranty	\$138.95
B***	PREMIUM TURBINE OVERHAUL : disassembly of turbine, cleaning, and replacement of all wear parts (bearings, o-rings, etc) using "euro-style" premium bearings; 8 month warranty	\$129.95
C	TURBINE REPLACEMENT : installation of complete turbine assembly:	
	• Bur tool style turbine, 6 month warranty	\$124.95
	• Pushbutton or lever-style turbine, 6 month warranty	\$159.95
	• "maintenance-free" Star 430 series, 6 month warranty	\$189.95
	• premium "maintenance-free" Star 430 series, 1 year warranty	\$249.95
D	OEM TURBINE REPLACEMENT : installation of complete original equipment manufacturer turbine assembly. Warranty per manufacturer	Varies by make and model

LOW SPEED HANDPIECES

OPT.	DESCRIPTION	
F	STANDARD MOTOR OVERHAUL : disassembly of motor components, cleaning, and replacement of all wear parts (bearings, o-rings, gears, etc); 6 month warranty	\$125.00 to 185.00
	Kavo motors start at	\$185.00
G	CONTRA-ANGLE/NOSECONE/HEAD REPAIRS : disassembly of handpiece, cleaning, and replacement of all wear parts; 6 month warranty	\$65.00 to 85.00
	Kavo contra-angles and nosecones start at	\$85.00
H	SCALER REPAIRS : disassembly of handpiece, cleaning, and replacement of all wear parts; 6 month warranty	\$125.00 to 185.00 by model
I	"Electric Attachments" repairs starting at estimate will be provided All manufacturers (Kavo, W&H, NSK/Brasseler, Bien Air, Star)	\$189.95

LOW SPEED TERMINOLOGY



HEAD

CONTRA-ANGLE

MOTOR

STRAIGHT NOSECONE

COUPLER

- E OTHER: choose this option for any repairs not listed above, including but not limited to:
- water/air line repairs
 - fiber optics replacement
 - handpiece reassembly
 - push button conversion
 - *all electric handpieces/motors, curing lights, Cavitrons, etc.*
- ESTIMATE REQUIRED - Please check box if you require an estimate or wish to be advised of recommended repair option

Estimate available upon request

ALL REPAIR OPTIONS INCLUDE A COMPREHENSIVE INTERNAL AND EXTERNAL CLEANING AND TESTING OF PERFORMANCE AND OUTPUT

*** PUSH BUTTON AND LEVER-STYLE TURBINES CONTAIN AN "AUTOCHUCKING" MECHANISM WHICH, IF ITS MECHANICAL INTEGRITY IS COMPROMISED, MAY FAIL TO PROPERLY GRIP THE CUTTING BUR SECURELY. DUE TO THE INTENSE SAFETY ISSUES INVOLVED, KAB Dental Inc WILL NOT REBUILD A TURBINE ON A FAULTY AUTOCHUCKING SPINDLE, AND WILL INSTEAD UPGRADE THE REPAIR TO OPTION "C", TURBINE REPLACEMENT